

PIL-BOB[®]



Makes it easy to get the pills out

Remove medication easily and hygienically without spillage

The Webstercare[®] Pil-Bob has been designed to make the removal of medication from Webster-pak[®] systems even easier.

The Pil-Bob is ideal for individuals who suffer from poor eyesight, arthritis or have poor dexterity, and suitable for aged care staff, self-medicating residents and carers helping people to administer medication. It has received favourable ergonomic review.

To use, insert the Pil-Bob tip into the blister compartment through the foil backing. Sweep it around the edge, allowing the medication to fall from the blister into the cup.

The advantages of Pil-Bob:

- Now manufactured with an antimicrobial additive proven to reduce the risk of bacterial cross contamination.
- Small and lightweight so it fits comfortably in your hand
- Durable plastic material makes it easy to clean and sterilise
- Centre tip with cup design prevents spilling or dropping medication
- Comes with OH&S instructions to ensure correct use (see overleaf)



Don't forget Cytotoxic Pil-Bob[®]

Used specifically to remove cytotoxic medications from Webster-pak[®] systems and help reduce the risk of cross contamination between medications and users.



"Innovative solutions to medication management"™



PIL-BOB®



Makes it easy to get the pills out

The Pil-Bob makes it even easier to remove pills hygienically and efficiently from Webstercare packs. **IMPORTANT:** Our experience shows if the Pil-Bob is used as indicated in the instructions it works very well. Please find below the instructions for use.

INSTRUCTIONS FOR USE:

STEP 1



Hold the Webster-pak folder in a horizontal (flat) position

STEP 2



Place the Pil-Bob behind the blister you wish to open

STEP 3



Insert the tip upwards through the foil backing

STEP 4



Sweep the blister with the serrated tip until all the pills drop out into the Pil-Bob cup

CAUTION: Swallowing pills directly from the Pil-Bob may result in injury due to the serrated tip. To avoid injury, transfer pills directly from the Pil-Bob into a disposable medication cup.

ORDER NOW from your Webster-pak pharmacy or direct from Webstercare®

TO ORDER OR FOR MORE INFORMATION ON THE PIL-BOB PLEASE CONTACT WEBSTERCARE

T (02) 9563 4900 | FREE CALL 1800 244 358
F (02) 9563 4955 | FREE FAX 1800 626 739
info@webstercare.com.au | www.webstercare.com.au

