

TECH TUESDAY

Webster-pak® 28

The Webstercare® Webster-pak 28 range has been designed to provide community pharmacies with a convenient, reliable and fully disposable medication management solution.

Featuring a durable cardboard frame for added strength and pack integrity, the Webster-pak 28 securely holds medications while offering the practicality of a disposable pack, making it an ideal option for customers who may not be able to return reusable folders to the pharmacy.

To suit the varying workflow and checking preferences of pharmacies, the range is available in three sealing options: Cold Seal, Clear Cold Seal and LowTemp™ foil. Each option offers its own unique benefits and is fully compatible with the MedsPro® Packing Robot, supporting efficient and streamlined medication packing processes.

Key features & benefits

Completely Disposable – Its disposable nature eliminates the need for returning any folders to the pharmacy – or having to chase them up from customers.

Pre-inserted blisters – The pre-inserted concertina blisters eliminate an additional assembly step when preparing the pack, helping to save valuable time during the packing process.

Pack can be hung on fridge and remain in sight – The ‘W’ at the top of the frame enables the pack to be hung on a fridge using the Webster-pak magnet, keeping medications well in sight, helping customers to remember to take each dose.

Burning questions


Q Which Webster-pak 28 should I choose?

A The Webster-pak 28 range is available in three different sealing options, each designed to suit different pharmacy workflows, packing preferences and checking requirements.

Webster-pak 28 Cold Seal Foil

Cold Seal is often selected when simplicity, and minimal equipment matter.


- Adhesive backing card with built in foil
- Most economic option
- No heat required - just peel and press to seal
- Faster packing workflow (no warming sealing device step)
- Less equipment dependency – useful for high-volume or mobile/remote packing



Webster-pak 28 LowTemp Foil

LowTemp is chosen when pharmacies want more flexibility and control over changes.


- Easiest to make adjustments or rework after sealing
- Works well where sealing devices are already part of the setup – you may already do heat seal packs



Webster-pak 28 Clear Seal Foil

Clear Seal is chosen mainly for accuracy and checking efficiency.

- Transparent backing allows medications to be seen from both sides
- Makes pharmacist checking faster and easier
- Helps reduce checking errors in busy environments
- Cold seal system (no heat required)



➤ [Click here to watch How to pack Webster-pak 28](#)

Have questions about our Webster-pak 28 range?

Get in touch with our customer service team on **1800 244 358** or email us at info@webstercare.com.au